

## GPR41 Antibody

Catalog No: #44332



Package Size: #44332-1 50ul #44332-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR41 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	GPR41 Antibody detects endogenous levels of total GPR41
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR41
Target Name	GPR41
Other Names	FFAR3, G-protein coupled receptor 41, FFA3R, Gpcr41, Gpr42, Free fatty acid receptor 3, G protein-coupled receptor 41, GPR41
Accession No.	Swiss-Prot#: O14843NCBI Gene ID:2865
Uniprot	O14843
GeneID	2865;
Calculated MW	38kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only