

Programmed cell death protein 2

Catalog No: #AP86600



Package Size: #AP86600-1 50ug #AP86600-2 100ug #AP86600-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Programmed cell death protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Zinc finger MYND domain-containing protein 7,Zinc finger protein Rp-8,RP8,ZMYND7
Accession No.	Uniprot ID:Q16342
Uniprot	Q16342
GeneID	5134;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDCD2

Note: This product is for in vitro research use only