

Polycomb group RING finger protein 6

Catalog No: #AP86602



Package Size: #AP86602-1 50ug #AP86602-2 100ug #AP86602-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Polycomb group RING finger protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mel18 and Bmi1-like RING finger,RING finger protein 134,MBLR,RNF134
Accession No.	Uniprot ID:Q9BYE7
Uniprot	Q9BYE7
GeneID	84108;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCGF6

Note: This product is for in vitro research use only