

Nuclear receptor subfamily 4 group A member 1

Catalog No: #AP86605



Package Size: #AP86605-1 50ug #AP86605-2 100ug #AP86605-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 4 group A member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Orphan nuclear receptor HMR, Orphan nuclear receptor NGFI-B, HMR
Accession No.	Uniprot ID:P51666
Uniprot	P51666
GeneID	403897;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NR4A1

Note: This product is for in vitro research use only