

NADH-ubiquinone oxidoreductase chain 5

Catalog No: #AP86606



Package Size: #AP86606-1 50ug #AP86606-2 100ug #AP86606-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NADH-ubiquinone oxidoreductase chain 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	NADH dehydrogenase subunit 5,MTND5,NADH5,ND5
Accession No.	Uniprot ID:Q9ZZ57
Uniprot	Q9ZZ57
GeneID	804484;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MT-ND5

Note: This product is for in vitro research use only