

Protein-serine O-palmitoleoyltransferase porcupine

Catalog No: #AP86607



Package Size: #AP86607-1 50ug #AP86607-2 100ug #AP86607-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein-serine O-palmitoleoyltransferase porcupine
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mPORC, Porc, Ppn
Accession No.	Uniprot ID:Q9JJJ7
Uniprot	Q9JJJ7
GeneID	53627;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Porcn

Note: This product is for in vitro research use only