

Peptidyl-prolyl cis-trans isomerase B

Catalog No: #AP86608



Package Size: #AP86608-1 50ug #AP86608-2 100ug #AP86608-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Peptidyl-prolyl cis-trans isomerase B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | PPlase B,CYP-S1,Cyclophilin B,Rotamase B,S-cyclophilin,SCYLP,CYPB |
| Accession No. | Uniprot ID:P23284 |
| Uniprot | P23284 |
| GeneID | 5479; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPIB

Note: This product is for in vitro research use only