

Protein phosphatase methylesterase 1

Catalog No: #AP86609



Package Size: #AP86609-1 50ug #AP86609-2 100ug #AP86609-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein phosphatase methylesterase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PME-1,PME1
Accession No.	Uniprot ID:Q58DN4
Uniprot	Q58DN4
GeneID	535390;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPME1

Note: This product is for in vitro research use only