

## PDZ and LIM domain protein 1

Catalog No: #AP86613



Package Size: #AP86613-1 50ug #AP86613-2 100ug #AP86613-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PDZ and LIM domain protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-terminal LIM domain protein 1,Elfin,LIM domain protein CLP-36,Clim1,Clp36
Accession No.	Uniprot ID:P52944
Uniprot	P52944
GeneID	54133;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pdim1

Note: This product is for in vitro research use only