

Oxidation resistance protein 1

Catalog No: #AP86619

Package Size: #AP86619-1 50ug #AP86619-2 100ug #AP86619-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxidation resistance protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Protein C7,C7
Accession No.	Uniprot ID:Q4V8B0
Uniprot	Q4V8B0
GeneID	117520;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Oxr1

Note: This product is for in vitro research use only