

Piezo-type mechanosensitive ion channel component 1

Catalog No: #AP86622

Package Size: #AP86622-1 50ug #AP86622-2 100ug #AP86622-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Piezo-type mechanosensitive ion channel component 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein FAM38A,Fam38a
Accession No.	Uniprot ID:E2JF22
Uniprot	E2JF22
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Piezo1

Note: This product is for in vitro research use only