

Paired immunoglobulin-like type 2 receptor beta

Catalog No: #AP86625



Package Size: #AP86625-1 50ug #AP86625-2 100ug #AP86625-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Paired immunoglobulin-like type 2 receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Activating receptor PILR-beta, Cell surface receptor FDFACT, Fdfact, Pilrb1
Accession No.	Uniprot ID:Q2YFS2
Uniprot	Q2YFS2
GeneID	170741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pilrb

Note: This product is for in vitro research use only