

## Secretory phospholipase A2 receptor

Catalog No: #AP86626



Package Size: #AP86626-1 50ug #AP86626-2 100ug #AP86626-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Secretory phospholipase A2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PLA2-R,180 kDa secretory phospholipase A2 receptor,M-type receptor
Accession No.	Uniprot ID:P49259
Uniprot	P49259
GeneID	317777;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PLA2R1

Note: This product is for in vitro research use only