

Serine/threonine-protein kinase PLK1

Catalog No: #AP86627



Package Size: #AP86627-1 50ug #AP86627-2 100ug #AP86627-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PLK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Polo-like kinase 1,PLK-1
Accession No.	Uniprot ID:Q2TA25
Uniprot	Q2TA25
GeneID	538238;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLK1

Note: This product is for in vitro research use only