

# Nucleotide-binding oligomerization domain-containing protein 1

Catalog No: #AP86629

Package Size: #AP86629-1 50ug #AP86629-2 100ug #AP86629-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nucleotide-binding oligomerization domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Caspase recruitment domain-containing protein 4,CARD4
Accession No.	Uniprot ID:Q9Y239
Uniprot	Q9Y239
GeneID	10392;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NOD1

Note: This product is for in vitro research use only