

Nucleolar GTP-binding protein 1

Catalog No: #AP86630



Package Size: #AP86630-1 50ug #AP86630-2 100ug #AP86630-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleolar GTP-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Chronic renal failure gene protein,GTP-binding protein NGB,Crfg,Nog1
Accession No.	Uniprot ID:Q99ME9
Uniprot	Q99ME9
GeneID	69237;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gtpbp4

Note: This product is for in vitro research use only