

Prostatic acid phosphatase

Catalog No: #AP86631



Package Size: #AP86631-1 50ug #AP86631-2 100ug #AP86631-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostatic acid phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5'-nucleotidase,5'-NT,Ecto-5'-nucleotidase,Fluoride-resistant acid phosphatase,FRAP,Thiamine monophosphatase,TMPase
Accession No.	Uniprot ID:P20646
Uniprot	P20646
GeneID	56780;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acpp

Note: This product is for in vitro research use only