

Peptidyl-prolyl cis-trans isomerase NIMA-interacting 1

Catalog No: #AP86632



Package Size: #AP86632-1 50ug #AP86632-2 100ug #AP86632-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase NIMA-interacting 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Peptidyl-prolyl cis-trans isomerase Pin1, PPIase Pin1, Rotamase Pin1
Accession No.	Uniprot ID:Q5BIN5
Uniprot	Q5BIN5
GeneID	535470;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PIN1

Note: This product is for in vitro research use only