

Pannexin-1

Catalog No: #AP86635



Package Size: #AP86635-1 50ug #AP86635-2 100ug #AP86635-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pannexin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MRS1,UNQ2529/PRO6028
Accession No.	Uniprot ID:Q96RD7
Uniprot	Q96RD7
GeneID	24145;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: PANX1

Note: This product is for in vitro research use only