

## Phosphoenolpyruvate carboxykinase, cytosolic [GTP]

Catalog No: #AP86636



Package Size: #AP86636-1 50ug #AP86636-2 100ug #AP86636-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphoenolpyruvate carboxykinase, cytosolic [GTP]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PEPCK-C,PPCK1
Accession No.	Uniprot ID:Q8HYZ4
Uniprot	Q8HYZ4
GeneID	282855;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PCK1

Note: This product is for in vitro research use only