

## Programmed cell death protein 4

Catalog No: #AP86638



Package Size: #AP86638-1 50ug #AP86638-2 100ug #AP86638-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Programmed cell death protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neoplastic transformation inhibitor protein,Nuclear antigen H731-like,Protein 197/15a,H731
Accession No.	Uniprot ID:Q53EL6
Uniprot	Q53EL6
GeneID	27250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PDCD4

Note: This product is for in vitro research use only