

Protein phosphatase 1 regulatory subunit 1A

Catalog No: #AP86640



Package Size: #AP86640-1 50ug #AP86640-2 100ug #AP86640-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein phosphatase 1 regulatory subunit 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Protein phosphatase inhibitor 1,I-1,IPP-1
Accession No.	Uniprot ID:Q8WMS3
Uniprot	Q8WMS3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PPP1R1A

Note: This product is for in vitro research use only