

# Phospholipase D2

Catalog No: #AP86645



Package Size: #AP86645-1 50ug #AP86645-2 100ug #AP86645-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospholipase D2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PLD 2,Choline phosphatase 2,PLD1C,Phosphatidylcholine-hydrolyzing phospholipase D2
Accession No.	Uniprot ID:O14939
Uniprot	O14939
GeneID	5338;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PLD2

Note: This product is for in vitro research use only