

Melanocyte protein PMEL

Catalog No: #AP86646



Package Size: #AP86646-1 50ug #AP86646-2 100ug #AP86646-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanocyte protein PMEL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	115 kDa melanosomal matrix protein, Melanocyte protein Pmel 17, Premelanosome protein, Silver locus protein homolog, MMP115, PMEL17, SILV
Accession No.	Uniprot ID:Q98917
Uniprot	Q98917
GeneID	396007;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PMEL

Note: This product is for in vitro research use only