

# Plasma cell-induced resident endoplasmic reticulum protein

Catalog No: #AP86648

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86648-1 50ug #AP86648-2 100ug #AP86648-3 1mg

## Description

Product Name	Plasma cell-induced resident endoplasmic reticulum protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Plasma cell-induced resident ER protein,pERp1,Proapoptotic caspase adapter protein
Accession No.	Uniprot ID:Q9D811
Uniprot	Q9D811
GeneID	69816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pacap

Note: This product is for in vitro research use only