

Prefoldin subunit 2

Catalog No: #AP86649



Package Size: #AP86649-1 50ug #AP86649-2 100ug #AP86649-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prefoldin subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PFD2,HSPC231
Accession No.	Uniprot ID:Q9UHV9
Uniprot	Q9UHV9
GeneID	5202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PFDN2

Note: This product is for in vitro research use only