

cGMP-specific 3',5'-cyclic phosphodiesterase

Catalog No: #AP86651



Package Size: #AP86651-1 50ug #AP86651-2 100ug #AP86651-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	cGMP-specific 3',5'-cyclic phosphodiesterase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	cGMP-binding cGMP-specific phosphodiesterase,CGB-PDE,PDE5
Accession No.	Uniprot ID:O76074
Uniprot	O76074
GeneID	8654;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDE5A

Note: This product is for in vitro research use only