

Platelet-derived growth factor subunit A

Catalog No: #AP86652



Package Size: #AP86652-1 50ug #AP86652-2 100ug #AP86652-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PDGF subunit A,PDGF-1,Platelet-derived growth factor A chain,Platelet-derived growth factor alpha polypeptide,Rpa1
Accession No.	Uniprot ID:P28576
Uniprot	P28576
GeneID	25266;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdgfa

Note: This product is for in vitro research use only