

Platelet-derived growth factor subunit B

Catalog No: #AP86653



Package Size: #AP86653-1 50ug #AP86653-2 100ug #AP86653-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PDGF subunit B,PDGF-2,Platelet-derived growth factor B chain,Platelet-derived growth factor beta polypeptide
Accession No.	Uniprot ID:Q6Q7I7
Uniprot	Q6Q7I7
GeneID	442986;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDGFB

Note: This product is for in vitro research use only