

[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial

Catalog No: #AP86654

Orders: order@signalwayantibody.com

Package Size: #AP86654-1 50ug #AP86654-2 100ug #AP86654-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pyruvate dehydrogenase kinase isoform 1,PDH kinase 1,PDHK1
Accession No.	Uniprot ID:Q15118
Uniprot	Q15118
GeneID	5163;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDK1

Note: This product is for in vitro research use only