

Period circadian protein homolog 3

Catalog No: #AP86655



Package Size: #AP86655-1 50ug #AP86655-2 100ug #AP86655-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Period circadian protein homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hPER3,Cell growth-inhibiting gene 13 protein,Circadian clock protein PERIOD 3,GIG13
Accession No.	Uniprot ID:P56645
Uniprot	P56645
GeneID	8863;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PER3

Note: This product is for in vitro research use only