

PCNA-associated factor

Catalog No: #AP86657



Package Size: #AP86657-1 50ug #AP86657-2 100ug #AP86657-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCNA-associated factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HCV NS5A-transactivated protein 9 homolog,p15PAF,Ns5atp9
Accession No.	Uniprot ID:Q6RIA2
Uniprot	Q6RIA2
GeneID	300795;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pclaf

Note: This product is for in vitro research use only