

Oncostatin-M-specific receptor subunit beta

Catalog No: #AP86659



Package Size: #AP86659-1 50ug #AP86659-2 100ug #AP86659-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oncostatin-M-specific receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Interleukin-31 receptor subunit beta,IL-31 receptor subunit beta,Osmrb
Accession No.	Uniprot ID:Q65Z14
Uniprot	Q65Z14
GeneID	310132;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Osmr

Note: This product is for in vitro research use only