

OTU domain-containing protein 4

Catalog No: #AP86660

Package Size: #AP86660-1 50ug #AP86660-2 100ug #AP86660-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OTU domain-containing protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HIV-1-induced protein HIN-1,HIN1,KIAA1046
Accession No.	Uniprot ID:Q01804
Uniprot	Q01804
GeneID	54726;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OTUD4

Note: This product is for in vitro research use only