

Orexin receptor type 2

Catalog No: #AP86662



Package Size: #AP86662-1 50ug #AP86662-2 100ug #AP86662-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Orexin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ox-2-R,Hypocretin receptor type 2
Accession No.	Uniprot ID:O43614
Uniprot	O43614
GeneID	3062;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HCRTR2

Note: This product is for in vitro research use only