

# Decorin

Catalog No: #AP86664



Package Size: #AP86664-1 50ug #AP86664-2 100ug #AP86664-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Decorin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Bone proteoglycan II,Dermatan sulfate proteoglycan-II,DSPG,PG-S2,PG40
Accession No.	Uniprot ID:Q01129
Uniprot	Q01129
GeneID	29139;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dcn

Note: This product is for in vitro research use only