

## Serine/threonine-protein kinase PAK 1

Catalog No: #AP86670



Package Size: #AP86670-1 50ug #AP86670-2 100ug #AP86670-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase PAK 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Alpha-PAK,p21-activated kinase 1,PAK-1,p65-PAK
Accession No.	Uniprot ID:Q08E52
Uniprot	Q08E52
GeneID	533729;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PAK1

Note: This product is for in vitro research use only