

Pyrraline-5-carboxylate reductase 2

Catalog No: #AP86673



Package Size: #AP86673-1 50ug #AP86673-2 100ug #AP86673-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyrraline-5-carboxylate reductase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	P5C reductase 2
Accession No.	Uniprot ID:Q17QJ7
Uniprot	Q17QJ7
GeneID	504987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PYCR2

Note: This product is for in vitro research use only