

## Aggrecan core protein

Catalog No: #AP86682



Package Size: #AP86682-1 50ug #AP86682-2 100ug #AP86682-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aggrecan core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cartilage-specific proteoglycan core protein,CSPCP,Chondroitin sulfate proteoglycan core protein 1,Chondroitin sulfate proteoglycan 1,AGC1,CSPG1,MSK16
Accession No.	Uniprot ID:P16112
Uniprot	P16112
GeneID	176;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ACAN

Note: This product is for in vitro research use only