

Phosphatidylinositol N-acetylglucosaminyltransferase subunit Y



Catalog No: #AP86683

Orders: order@signalwayantibody.com

Package Size: #AP86683-1 50ug #AP86683-2 100ug #AP86683-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Phosphatidylinositol N-acetylglucosaminyltransferase subunit Y
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Phosphatidylinositol-glycan biosynthesis class Y protein,PIG-Y
Accession No.	Uniprot ID:P0C1N9
Uniprot	P0C1N9
GeneID	514374;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PIGY

Note: This product is for in vitro research use only