

## Pituitary homeobox 2

Catalog No: #AP86684



Package Size: #AP86684-1 50ug #AP86684-2 100ug #AP86684-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pituitary homeobox 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ALL1-responsive protein ARP1,Homeobox protein PITX2,Paired-like homeodomain transcription factor 2,RIEG bicoid-related homeobox transcription factor,Solurshin,ARP1,RGS,RIEG,RIEG1
Accession No.	Uniprot ID:Q99697
Uniprot	Q99697
GeneID	5308;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PITX2

Note: This product is for in vitro research use only