

ProSAAS

Catalog No: #AP86686



Package Size: #AP86686-1 50ug #AP86686-2 100ug #AP86686-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ProSAAS
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IA-4,Proprotein convertase subtilisin/kexin type 1 inhibitor,Proprotein convertase 1 inhibitor,pro-SAAS
Accession No.	Uniprot ID:Q9QXV0
Uniprot	Q9QXV0
GeneID	30052;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pcsk1n

Note: This product is for in vitro research use only