

## Protein disulfide-isomerase

Catalog No: #AP86687



Package Size: #AP86687-1 50ug #AP86687-2 100ug #AP86687-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein disulfide-isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDI, Cellular thyroid hormone-binding protein, Endoplasmic reticulum resident protein 59, ER protein 59, Prolyl 4-hydroxylase subunit beta, p55, Pdia1
Accession No.	Uniprot ID:P09103
Uniprot	P09103
GeneID	18453;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: P4hb

Note: This product is for in vitro research use only