

Period circadian protein homolog 3

Catalog No: #AP86688



Package Size: #AP86688-1 50ug #AP86688-2 100ug #AP86688-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Period circadian protein homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rPER3,Circadian clock protein PERIOD 3
Accession No.	Uniprot ID:Q8CJE2
Uniprot	Q8CJE2
GeneID	78962;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Per3

Note: This product is for in vitro research use only