

Aggrecan core protein

Catalog No: #AP86693



Package Size: #AP86693-1 50ug #AP86693-2 100ug #AP86693-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aggrecan core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cartilage-specific proteoglycan core protein,CSPCP,Agc,Agc1
Accession No.	Uniprot ID:Q61282
Uniprot	Q61282
GeneID	11595;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Acan

Note: This product is for in vitro research use only