

Plasma glutamate carboxypeptidase

Catalog No: #AP86695



Package Size: #AP86695-1 50ug #AP86695-2 100ug #AP86695-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasma glutamate carboxypeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hematopoietic lineage switch 2 related protein,Hls2-rp,Liver annexin-like protein 1,LAL-1
Accession No.	Uniprot ID:Q6IRK9
Uniprot	Q6IRK9
GeneID	58952;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pgcp

Note: This product is for in vitro research use only