

Prolactin-inducible protein homolog

Catalog No: #AP86697

Package Size: #AP86697-1 50ug #AP86697-2 100ug #AP86697-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolactin-inducible protein homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	14 kDa submandibular gland protein, SMGP, Gross cystic disease fluid protein 15, GCDFP-15, Prolactin-induced protein
Accession No.	Uniprot ID:P02816
Uniprot	P02816
GeneID	18716;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pip

Note: This product is for in vitro research use only