

Poly(A) RNA polymerase, mitochondrial

Catalog No: #AP86701



Package Size: #AP86701-1 50ug #AP86701-2 100ug #AP86701-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Poly(A) RNA polymerase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PAP,PAP-associated domain-containing protein 1,Polynucleotide adenylyltransferase,Terminal uridylyltransferase 1,TUTase 1,mtPAP,PAPD1
Accession No.	Uniprot ID:Q9NVV4
Uniprot	Q9NVV4
GeneID	55149;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MTPAP

Note: This product is for in vitro research use only