

Putative phospholipase B-like 2

Catalog No: #AP86704



Package Size: #AP86704-1 50ug #AP86704-2 100ug #AP86704-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Putative phospholipase B-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	76 kDa protein,p76,LAMA-like protein 2,Lamina ancestor homolog 2,Phospholipase B domain-containing protein 2
Accession No.	Uniprot ID:Q8NHP8
Uniprot	Q8NHP8
GeneID	196463;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLBD2

Note: This product is for in vitro research use only