

Plastin-3

Catalog No: #AP86706



Package Size: #AP86706-1 50ug #AP86706-2 100ug #AP86706-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Plastin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	T-plastin
Accession No.	Uniprot ID:P13797
Uniprot	P13797
GeneID	5358;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLS3

Note: This product is for in vitro research use only